



Medley Ergo: cost effective patient care

The **Medley Ergo** range from Invacare offers the perfect balance between budget, functionality and quality requirements. The four-sectioned profiling mattress support allows for a range of height combinations, providing an ideal position for the carer and comfort for the client. The newly optimised mattress platform sections reduce shear and friction when the bed is being adjusted and prevents the client from sliding downwards.

The **Medley Ergo Low** offers the same benefits as the **Medley Ergo** but lowers to a 210mm height to maximise safety for vulnerable clients who are prone to unexpected bed leaving.



Trusted Solutions, Passionate People





Optimised dimensions of the mattress platform

The mattress platform sections boast optimal dimensions designed in accordance with anthropometric statistics and pressure mapping recommendations. The **Medley Ergo** range offers:

- Increased comfort when sitting
- Reduced shear and friction due to sliding downwards during profiling
- Also suitable for shorter clients

Improved functionality and efficiency

The *Invacare Medley Ergo* range has been specially designed for quick transportation to meet the needs of the transient working environments of Community Care provision. The **Medley Ergo** can be supplied on a transporter kit so the bed can easily be moved and installed. In addition, any assembling can be done without the use of tools or technicians.

Height is no longer an issue

The flexible design of the **Medley Ergo** range allows the bed to be mounted at a low height, or a higher position if preferred, creating a range of working heights for the carer and

transfer heights for the client. This is facilitated by the unique feature of the bed end, which enables the mattress support to be mounted in both upper and lower positions.

Features and Options



Full length side rails

Available with wooden or steel full length side rails. Wooden extendable side rails (150mm) are also available.



Scala 2 range

A range of collapsible three quarter length steel side rails available in various heights to cover standard to deep mattress heights.



Safety in and out of bed

In the event of a client falling from the bed, a foldable mat can be placed on the floor and attached to the **Medley Ergo Low** for added safety. Order code: 1494608.

Options



Medley Ergo Low (21cm lowest height)

Medley Ergo Low

The height of the **Medley Ergo Low** bed offers increased security for clients enabling them to sleep with real peace of mind. The **Medley Ergo Low** sits low to the floor, (210mm in height) which, combined with a floor mattress, reduces harm from falls. Carers also have the option to raise the bed to a more suitable height when necessary.

The **Medley Ergo Low** makes it easier for clients to get in and out of bed, providing more independence for those who can self transfer.

The **Medley Ergo Low** is dual height adjustable, offering a 210-610mm and 280-680mm height option to suit a wider range of needs.



Medley Ergo Select

Homely design

The **Medley Ergo Select** offers an attractive, homely design. Also available as a retro-fit to the **Medley Ergo** and **Medley Ergo Low**, this aesthetic bed end cover fits with precision onto the bed actuators and can be removed at any time.

The leg supports are equipped with a Rastofix to enable individual positioning for the legs. **Medley Ergo Select** can easily be dismantled and stored, using a transport adapter for moving and storage.

Features and options



Easy to transport

Medley Ergo beds are easy to assemble and dismantle and can be transported using transporter kits (optional).



Lifting Pole

Can be adjusted in both height and depth (optional).



Select design for Medley Ergo

The homely design of **Medley Ergo Select** can also be retro-fitted (optional).



Double height positions for bed ends

Adjustable height range for **Medley Ergo** of 330-730mm or 400-800mm.

Medley Ergo Low 210-610mm or 280-680mm.



Optimised dimensions of mattress platform

According to anthropometrical measurements, the optimised dimensions meets the body size of a wider number of clients and reduces shear and friction.

Features and Options



Mattress support extension
Extension kit for increasing mattress support length by 150mm.



Transport fitting
For easy storage and transport.



Hand controls
Available with or without lockable functions.

Technical data

						Complies to the new safety standard, IEC 60601-2-52
Medley Ergo	Width outside: 1020 mm Width inside: 900 mm	Length outside: 2240 mm Length inside: 2080 mm	Height adjustment: 400 - 800 mm / 330 - 730 mm	Mattress support dimensions: 800-230-300-670 mm	0 - 70° Backrest angle	
Medley Ergo Low	Width outside: 1020 mm Width inside: 900 mm	Length outside: 2240 mm Length inside: 2080 mm	Height adjustment: 280 - 680 mm / 210 - 610 mm	Mattress support dimensions: 800-230-300-670 mm	0 - 70° Backrest angle	0 - 24° Thigh angle

Complies to the new safety standard, IEC 60601-2-52

For mattress compatibility, please refer to the user manual which can be downloaded from www.invacare.co.uk.

					Invacare Limited Pencoed Technology Park - Pencoed - Bridgend CF35 5AQ - United Kingdom Tel: +44 1 656 776 222 Fax: +44 1656 776 220 E-mail: uk@invacare.com Sales Order E-mail: ordersuk@invacare.com www.invacare.co.uk
Medley Ergo	0 - 11° Legrest lift	11° Tilt option ¹	64 kg	SWL 180 kg (28.5 stone) Max. user weight 145 kg (23 stone)	
Medley Ergo Low	0 - 11° Legrest lift	11° Tilt option ¹	64 kg	SWL 180 kg (28.5 stone) Max. user weight 145 kg (23 stone)	Heaviest part of the product: 17 kg

^{1/} Please note: the tilt function is only used for clients under medical supervision. If there is a risk of clients using the tilt function without supervision, please use a lockable hand control.

NB/ The bed must not be used by patients under 12 years of age, or by patients with body size equivalent to an average 12 year old or smaller.

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS -
County Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

Frame colour and wood work



Light Grey



Beech

Warranty

3 years frame
3 years actuator
2 years handset

© 2016 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Medley Ergo - UK - 03/2016 - 1554900

The *Invacare Medley Ergo* is part of the Safe Patient Handling range.



safe patient handling



Yes, you can.®